




Investigations Into Brownian Motion:

Optical Trapping

MIT Department of Physics
10/22/21


Vedang Lad

Partner: Bryan Sperry



What is Optical Trapping

- Piconewton sized forces on micron sized dielectric objects with infrared laser
- Revolutionary in Biophysics
- Pick up, move around, and interact with small biological systems



Arthur Ashkin

Goal: Determine the Boltzmann Constant

How has it been done before?



- Dielectric-constant gas thermometry
 - Gas response to electric fields
- Acoustic thermometry
 - Speed of sound in a gas is directly dependent on its temperature

How will we do it?



- Brownian motion and power spectral density using optical trapping

The Physics of Optical Trapping

- Momentum
- Scattered light
creates a scattering
force
- Refracted light
creates gradient
force



Set Up



Trap in X-Y
Focus laser with Z

Math of Optical Trapping

β hydrodynamic drag coefficient

α harmonic optical trap of stiffness

d diameter of bead

η viscosity of the medium

Wiener-Khinchin theorem to define a
“power spectral distribution” function

Equipartition theorem...

$$\beta x' + \alpha x = F(t)$$

$$\beta = 3\pi\eta d$$

$$S_{xx}(f) = \frac{K_B T}{\pi^2 \beta (f^2 + f_0^2)}$$

$$f_0 = \frac{\alpha}{2\pi\beta}$$

Step One: Obtain Calibration






- Problem?
 - QPD - millivolts
 - Photodiode
 - Piezo - volts
 - Strain gauge
- How do we extract the position?
 - Pass through over fixed bead sample!

Part 1: Relationship between QPD and Piezos

Reducing systematic errors:

- Three trials for each X and Y at each laser current for the fixed bead sample
- Laser intervals for experiment: 100, 125, 150, 175, 200








Part 2: Relationship between Piezo and Position



- Four different points
- Preset at two



The Complete Conversion from QPD to Position



$$QPD(mV) \times \frac{Piezo (V)}{QPD(mV)} \times \frac{Position (\mu m)}{Piezo (V)} = \mu m$$

$[Part\ 1: Slope]^{-1}$
 $Part\ 2: Slope$

Step Two: Power Spectral Distribution

- Fourier transform QPD data using signal processing trick
- Fit
 - Extract boltzmann constant
 - Extract stiffness
- Function to fit:


$$S_{xx}(f) = \frac{A}{f^2 + f_0^2}$$

$$\beta = 3\pi\eta d$$



$$S_{xx}(f) = \frac{K_B T}{\pi^2 \beta (f^2 + f_0^2)}$$

$$f_0 = \frac{\alpha}{2\pi\beta}$$

Power Spectral Distribution Plot



Fitted Power Spectrum Density: QPD X and QPD Y




Putting It All Together


The diagram illustrates the process of substituting known parameters into a general equation for the spectral density $S_{xx}(f)$. On the left, the general equation is $S_{xx}(f) = \frac{A}{f^2 + f_0^2}$. Three teal arrows originate from the right side of this equation, pointing to the definitions of the parameters A , f_0 , and β respectively. The top arrow points to $\beta = 3\pi\eta d$. The middle arrow points to $S_{xx}(f) = \frac{K_B T}{\pi^2 \beta (f^2 + f_0^2)}$, where the entire right-hand side of this equation is circled in teal. The bottom arrow points to $f_0 = \frac{\alpha}{2\pi\beta}$.

$$S_{xx}(f) = \frac{A}{f^2 + f_0^2}$$
$$\beta = 3\pi\eta d$$
$$S_{xx}(f) = \frac{K_B T}{\pi^2 \beta (f^2 + f_0^2)}$$
$$f_0 = \frac{\alpha}{2\pi\beta}$$

Trap Stiffness for X and Y



Determination of Boltzmann Constant for X and Y



Grand Finale: Boltzmann Constant

$$k_B = QPD \text{ X: } 1.79 \times 10^{-23} \pm 1.10 \cdot 10^{-23} \frac{J}{K}$$

$$k_B = QPD \text{ Y: } 1.72 \times 10^{-23} \pm 3.94 \cdot 10^{-24} \frac{J}{K}$$

Error Discussion

- Heating due to laser $T = 293 \pm 5 \text{ K}$ in reality. ● 10%
- Piezo electronics (strain gauge) ● 5%
- Lens and coverslip ● 3%
- Other ● 1%
- Propagation (statistical)

Conclusion

$$k_B = QPD\ X: 1.79 \times 10^{-23} \pm 1.10 \cdot 10^{-23} \text{ (stat)} \pm 0.34 \times 10^{-23} \text{ (sys)} \frac{J}{K}$$

$$k_B = QPD\ Y: 1.72 \times 10^{-23} \pm 0.39 \cdot 10^{-23} \text{ (stat)} \pm 0.32 \times 10^{-23} \text{ (sys)} \frac{J}{K}$$


- Calibration but inherent dependence on electronics
- Measure active temperature of water
- No way to quantify interactions between stage, slip, oil, and lens

They were within error of each other different from the true value of $1.38 \times 10^{-23} \frac{J}{K}$

Resources

- [1] J. L. Staff, Photoelectric effect lab guide (2018), jLab E-Library, URL <http://web.mit.edu/8.13/www/JLExperiments/JLExp51.pdf>
- [2] UC. B. Staff, Advanced Experimentation Laboratory (2018), University of California, Berkeley, URL <http://experimentationlab.berkeley.edu/sites/default/files/writeups/OTZ.pdf>
- [3] P. Bevington and D. Robinson, Data Reduction and Error Analysis for the Physical Sciences (McGraw-Hill, 2003).

**Thank
You**



Average Chi Square Probability:
0.702 , 0.771

Standard Deviation: 0.261, 0.229

